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Thriller® Thermoshaker Incubator / Thermo-Inkubationsmischer

CONTENT

WARRANTY	2
PACKAGING LIST	2
INSTALLATION	2
OPERATING ELEMENTS	3
OPERATING THE THRILLER®	4
Direct operation program P0	4
Programs P1 to P9	5
Use thermostat without mixing activities	5
Use thermostat with mixing activities	5
Mixing without incubating facility	6
Mixing at intervals	6
Short mix	6
Storing parameters	6
CLEANING, MAINTENANCE AND REPAIR	7
TECHNICAL DATA	7
TECHNICAL SUPPORT AND ORDERING INFORMATION	8
DISPOSAL INSTRUCTION	9

WARRANTY

PEQLAB guarantees that the system you have received has been thoroughly tested and meets its published specification.

However, immediately upon arrival, please check carefully that the shipment is complete and has not been damaged in transit. For missing parts or to report any kind of damage, please contact PEQLAB (see 'TECHNICAL SUPPORT AND ORDERING INFORMATION'). Please retain all packaging materials until the delivery has been completely checked since this will speed up the return of goods if required and reduce environmental impact. Any form of returns, replacements or credit notes must be agreed in advance by PEQLAB.

For the Thriller® Thermoshaker Incubator and the thermoblocks PEQLAB guarantees a warranty period of 24 months if the products have been used solely according to the instruction manual and if not agreed differently. No liability is accepted for loss or damage arising from incorrect use. PEQLAB's liability is limited to the repair or replacement of the unit or refund of the purchase price, at PEQLAB's discretion. PEQLAB is not liable for any consequential damages.

PEQLAB reserves the right to alter the technical specifications without prior notice. This will enable us to implement developments as soon as they arise.

PACKAGING LIST

Thriller® Thermoshaker Incubator without block
1 main power cable
1 hex driver
1 instruction manual

INSTALLATION

The Thriller® should be placed on a firm support plate. Sufficient ventilation must be ensured for perfect functional operation of the Thriller®. For this reason, the Thriller® should always have a few centimeters of space between the unit and neighboring instruments or the wall.

Remove the packaging and place the Thriller® on a firm support plate. Remove the packaging of the block. The block can be fitted to the Thriller® in one direction only. The orientation is defined by the measuring probe stylus.

Place the block on the Thriller® in the correct orientation and use the hex driver enclosed to tighten the hexagon head cap screws.

Use the main cable supplied to connect the Thriller® to the main supply. The connector and the main switch are located on the rear side of the Thriller®.

OPERATING ELEMENTS


Figure 1: Operating elements of the Thriller®

Start/Stop button	To start or stop the program displayed
Prog button	To change from P0 to the stored programs P1 to P9
Save button	To save the programs P1 to P9
Short mix button	The Thriller® will continue to mix at the intensity shown as long as the button is pressed
+/- buttons	To increase/decrease the values
←/→	To change to other parameters

Programs

The Thriller® has a direct operation program P0. The parameters will be saved automatically after input. If the setpoint temperature is changed, the instrument will start heating or cooling immediately without having to start the program. If time or mixing functions are to be considered, the program must be activated by pressing the Start/Stop button.

Another nine program sites (P1 to P9) with nine sub-steps each (S1 to S9) are available to save programs. The Start/Stop button must be pressed whenever a program is to be activated. Changes are saved by pressing the Save button; they are not saved automatically.

The Thriller® could be programmed "infinite", when total time is set to the maximum. The maximum set point is 99h 59min, after this you can choose "inf" for infinite.

OPERATING THE THRILLER®

Switch on the Thriller® at the rear side of the instrument, the display will show the values of the direct operation program P0 set last.

Set. temp./Act. temp.	Total time	Mixing intensity
60/22	01h00m	450 rpm
00:10m	00:05m	P0
Mixing activity	Mixing break	Prog No./Prog step

Figure 2: Display of the Thriller®

Direct operation program P0

- After the instrument has been switched on the cursor will start flashing at the setpoint temperature.
- Use the + or – button to set the requested temperature. The Thriller® will start the heating or cooling process automatically.
- Use the → key to skip to the total time. Use the + or – button to set the total time. The Thriller® could be programmed “infinite”, when total time is set to the maximum. The maximum set point is 99h 59min, after this you can choose “inf” for infinite.
- Use the → key to skip to the mixing intensity. Use the + or – button to set the mixing intensity.
- Use the → key to skip to the mixing time. Use the interval function to enter a continuous change of mixing activity and mixing break. If no interval is required, set these parameters to OFF. If an interval function is required, use the + or – button to set the time of the mixing activity. Use the → button to skip to the mixing break, and then use the + or – button to set the time of the break. If the mixing time has been set to OFF, the mixing intensity set is used for continuous mixing.

Programs P1 to P9

- Press the Prog button to skip to the programs saved or use one of the arrow keys to move the cursor to the program number.
- Use the + or – button to select the requested program. Always start at program step S1. Later program steps will be taken into consideration only, if the previous step contains a programmed time. Not equal to zero or infinite. Subsequently the next program step will be started automatically.
- Use the arrow keys to move the cursor to the setpoint temperature.
- Use the + or – button to set the requested temperature.
- Use the → key to skip to the total time. Use the + or – button to set the total time. The Thriller® could be programmed “infinite”, when total time is set to the maximum. The maximum set point is 99h 59min, after this you can choose “inf” for infinite.
- Use the → key to skip to the mixing intensity. Use the + or – button to set the mixing intensity.
- Use the → key to skip to the mixing time. Use the interval function to enter a continuous change of mixing activity and mixing break. If no interval is required, set these parameters to OFF. If an interval function is required, use the + or – button to set the duration of the mixing activity.
- Use the → button to skip to the mixing break, and then use the + or – button to set the duration of the break. If the mixing time has been set to OFF, the mixing intensity set is used for continuous mixing.
- Save any program changes by pressing the Save button.
- Start the program by pressing the Start/Stop button.

In case of programs with several steps, move the cursor to the step number, and use the + button to skip to S2. Enter parameters as described above. The program ends at the last program step at which the total time has been set.

Use thermostat without mixing activities

Set the temperature and the total time; set the mixing intensity to 0. Start the program. The Thriller® will temperature until the time has run down. A signal is sounded as soon as the time has run down.

Use thermostat with mixing activities

Set the temperature, the total time and the mixing intensity. Start the program. The Thriller® will temperature and mix until the time has run down. A signal is sounded as soon as the time has run down.

Mixing without incubating facility

Set the temperature to OFF. The Thriller® will cool down or heat up until the block has reached 25 °C. Subsequently the block will adopt the ambient temperature.

Set the total time and the mixing intensity and start the program. A signal is sounded as soon as the time has run down.

Mixing at intervals

Set the setpoint temperature or set the temperature to OFF. Set the total time and the mixing intensity. Set the requested mixing activity (e.g. 0:10m = 10 seconds). Subsequently set the time of the mixing break (e.g. 0:05m = 5 seconds). Start the program. The Thriller® will start operating at a change of mixing activity (0:10m) and mixing break (0:05m) in succession. After the total time has run down, the program will stop.

Short mix

The Thriller® mixes at the intensity shown in the display as long as the button is pressed.

Storing parameters

Changes in the programs P1 to P9 are saved by pressing the Save button. If a program is quit without saving, the query will appear whether the changes shall be saved ('save changes?'). The changes are saved by pressing the + button, whereas pressing the – button will cause the program to quit without saving.

Changes in the direct operation program P0 will be saved automatically.

CLEANING, MAINTENANCE AND REPAIR

Please make sure to disconnect the instrument from the main power if cleaning is desired. We recommend the use of a mild soap solution to clean the housing. If sample liquid has run into the boreholes of the blocks, they should be cleaned with a mild soap solution immediately and rinsed with distilled water. Mild alcohol-containing disinfectants may be used for decontamination. Maintenance of the Thriller® is not required.

In case a repair gets necessary the device must be sent to the manufacturer. **The surface of the device should be decontaminated with a lint free cloth soaked with 70 % ethanol before sending.** The drill holes of the thermoblock shall be cleaned with a 1.5 % hypochlorite solution with the help of a cotton bud. **Please fill in and sign the decontamination certificate and send it to the manufacturer together with the device.**

TECHNICAL DATA

Temperature control range:	4 °C (at 20 °C RT) to 99 °C in 1 °C steps
Temperature accuracy:	± 0.5 °C

Heating rate:	approx. 5 °C/min
Cooling rate:	4 °C/min 99 °C to 20 °C 1 °C/min 20 °C to 4 °C

The heating and cooling rates may differ when using fluid-filled reaction tubes or at higher ambient temperature.

Mix intensity in steps of 50 rpm:	
0.5 ml thermoblock	300 to 1500 rpm
1.5 and 2.0 ml thermoblock	300 to 1400 rpm
0.2 ml thermoblock	300 to 1400 rpm
Microplate thermoblock	300 to 1400 rpm

Mixing stroke:	3 mm
Programs:	1 direct operation program (P0) 9 programs with 9 steps each (P1 – P9)

Max. total time:	99 h 59 m or infinite
Max. mixing activity/pause (interval function):	99 m 59 s

Thermoblocks:	24 x 1.5 ml block 24 x 0.5 ml block 24 x 2.0 ml block 96 x 0.2 ml block Microplate block
---------------	--

Power supply:	115/230 V, 50 - 60 Hz
Power requirement:	125 W
Computer link:	RS232
Size:	19 x 19 x 32.5 cm
Weight:	5.9 kg incl. thermoblock
Max. humidity:	70 %
Ambient temperature:	4 - 35 °C

TECHNICAL SUPPORT AND ORDERING INFORMATION

For technical questions and more detailed information on PEQLAB's products please visit www.peqlab.com to find the respective contact person.

Please have the unit's serial number (located on the bottom panel of the instrument) available when calling. Should an item require return to PEQLAB for service, a decontamination form must be completed first by the user. Items sent without decontamination form will not be accepted.

All returns must be pre-approved by PEQLAB!

Thriller® Thermoshaker Incubator	91-7010
Thriller® Thermoshaker Incubator, 110 V version for USA	91-7010-US

Thermoblock 24 x 0.5 ml micro test tubes	91-7110
Thermoblock 24 x 1.5 ml micro test tubes	91-7111
Thermoblock 24 x 2.0 ml micro test tubes	91-7112
Thermoblock 96 x 0.2 ml micro test tubes or plates	91-7113
Thermoblock Microplate (max. suited plate height = 18 mm)	91-7114

MplTI SafeSeal Tubes 0.65 ml	
natural, 1 bag of 1000	81-11150
natural, Low Binding, 2 bags of 500	81-11330

MplTI SafeSeal Tubes 1.7 ml	
natural, 1 bag of 1000	81-11510
natural, Low Binding, 2 bags of 500	81-11720


MplTI SafeSeal Tubes 2.0 ml	
natural, 1 bag of 1000	81-12000
natural, Low Binding, 2 bags of 500	81-12180

96 well PCR Plates	
Standard 96 well, non-skirted, 25 plates	82-0600-A
Low Profile 96 well, non-skirted, 25 plates	82-0700-A
Semi-skirted 96 well, 25 plates	82-0900-A
QPCR semi-skirted 96 well, 25 plates	82-1400-A
Skirted 96 well, 25 plates	82-0800-A

DISPOSAL INSTRUCTION

This product is subject to the WEEE disposal regulations and cannot be disposed of with regular waste.

Please contact PEQLAB Biotechnologie GmbH or your local distributor for disposal of Thriller® Thermoshaker Incubator. Waste electronic equipment is not useless waste. Environmentally friendly disposal can help to retrieve valuable resources. Thereby you contribute to the protection of the environment and human health.

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